

572

(FILE 'HOME' ENTERED AT 09:01:00 ON 30 SEP 2005)

INDEX 'ADISCTI, ADISINSIGHT, ADISNEWS, AGRICOLA, ANABSTR, ANTE, AQUALINE, AQUASCI, BIOBUSINESS, BIOCOMMERCE, BIOENG, BIOSIS, BIOTECHABS, BIOTECHDS, BIOTECHNO, CABA, CANCERLIT, CAPLUS, CEABA-VTB, CEN, CIN, CONFSCI, CROPB, CROPU, DDFB, DDFU, DGENE, DISSABS, ...' ENTERED AT 09:01:10 ON 30 SEP 2005
SEA (FLUORESCENT (3N) DYE) (P) (DSRNA OR (DOUBLE (A) STRANDED

0* FILE ADISNEWS
0* FILE ANTE
0* FILE AQUALINE
2 FILE AQUASCI
0* FILE BIOCOMMERCE
1* FILE BIOENG
7 FILE BIOSIS
51* FILE BIOTECHABS
51* FILE BIOTECHDS
5* FILE BIOTECHNO
17 FILE CAPLUS
1* FILE CEABA-VTB
0* FILE CIN
1 FILE DDFU
66 FILE DGENE
1 FILE DISSABS
1 FILE DRUGU
1 FILE EMBASE
5* FILE ESBIODASE
3* FILE FEDRIP
0* FILE FOMAD
0* FILE FOREGE
0* FILE FROSTI
0* FILE FSTA
50 FILE IFIPAT
0* FILE KOSMET
2 FILE LIFESCI
5 FILE MEDLINE
0* FILE NTIS
0* FILE NUTRACEUT
5* FILE PASCAL
0* FILE PHARMAML
2 FILE PROMT
2 FILE SCISEARCH
3 FILE TOXCENTER
184 FILE USPATFULL
13 FILE USPAT2
0* FILE WATER
34 FILE WPIDS
34 FILE WPINDEX

L1 QUE (FLUORESCENT (3N) DYE) (P) (DSRNA OR (DOUBLE (A) STRANDED

FILE 'MEDLINE, CAPLUS, BIOSIS, BIOTECHNO, SCISEARCH' ENTERED AT 09:07:24
ON 30 SEP 2005

L2 36 S L1
L3 26 DUP REM L2 (10 DUPLICATES REMOVED)
L4 15 S L3 AND PY<2002
L5 1 S (FLUORESCENT (3N) DYE) (S) (BIND### OR ATTACH#####) (S) (DS
L6 7 S (FLUORESCENT (3N) DYE) (S) (BIND### OR ATTACH#####) (S) (DS
L7 8 S ((CYANIN! OR FLUORESCENT) (3N) DYE) (S) (BIND### OR ATTACH#
L8 8 S ((CYANIN? OR FLUORESCENT) (3N) DYE) (S) (BIND### OR ATTACH#
L9 7 DUP REM L8 (1 DUPLICATE REMOVED)
L10 5 S (CYANIN? (3N) DYE) (S) (DSRNA OR (DOUBLE (A) STRANDED OR DO
L11 5 DUP REM L10 (0 DUPLICATES REMOVED)
L12 4 S L11 NOT L9

L4 ANSWER.7 OF 15 CAPLUS COPYRIGHT 2005 ACS on STN
 AN 1998:501215 CAPLUS
 DN 129:145614
 TI Fluorescent intercalative dye assay method for nucleic acid sequences
 IN Ishiguro, Takahiko; Saitoh, Juichi
 PA Tosoh Corp., Japan
 SO Eur. Pat. Appl., 14 pp.
 CODEN: EPXXDW
 DT Patent
 LA English
 FAN.CNT 2

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
PI	EP 855447	A2	19980729	EP 1998-300481	19980123 <--
	EP 855447	A3	19990414		
	EP 855447	B1	20040331		
	R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO				
	JP 10201476	A2	19980804	JP 1997-10996	19970124 <--
	US 6063572	A	20000516	US 1998-12573	19980123 <--
	EP 1400598	A1	20040324	EP 2003-27996	19980123
	R: DE, FR, GB, IT				
PRAI	JP 1997-10996	A	19970124		
	EP 1998-300481	A3	19980123		

AB A method of assay of a specific nucleic acid anticipated in a sample, which comprises: (1) a DNA-producing step which involves production of a **double-stranded** DNA having a promoter sequence for an **RNA** polymerase and the nucleic sequence of the specific nucleic acid (the specific nucleic acid sequence) downstream from the promoter sequence by using the specific nucleic acid in the sample as a template, and (2) an **RNA**-producing and measuring step which involves production of a single-stranded **RNA** having the specific nucleic acid sequence by the **RNA** polymerase and measurement of the single-stranded **RNA**, and (3) wherein the **RNA**-producing and measuring step is initiated by adding at least the **RNA** polymerase, ribonucleoside triphosphates and a probe which is labeled with a **fluorescent intercalative dye** and is complementary to the single-stranded **RNA** to the reaction solution after the DNA producing step, involves measurement of the fluorescence intensity of the reaction solution, and is carried out at a constant temperature, and does not involve denaturing and annealing procedure for hybridization or separation of the probe which has not hybridized with the single-stranded **RNA** produced. The method is exemplified by anal. of recombinant RNA from human cytomegalovirus. The promoter primer used is designed to have at least the promoter sequence for RNA polymerase. A **fluorescent intercalative dye** such as YO-271 is linked to the oligonucleotide by a covalent bond.

L4 ANSWER 8 OF 15 CAPLUS COPYRIGHT 2005 ACS on STN
 AN 1992:102021 CAPLUS
 DN 116:102021
 TI Flow cytometric methods for RNA content analysis
 AU Darzynkiewicz, Zbigniew
 CS Cancer Res. Inst., New York Med. Coll., Valhalla, NY, 10595, USA
 SO Methods (San Diego, CA, United States) (1991), 2(3), 200-6
 CODEN: MTHDE9; ISSN: 1046-2023
 DT Journal
 LA English
 AB Anal. of the RNA content of individual cells provides information on their translational capacity that varies with cell differentiation and proliferation. In fact, because the rates of cell growth and proliferation are coupled, cellular or nuclear (nucleolar) RNA content indirectly serves as a marker of cell proliferation. The most common use of RNA measurement is in the discrimination of cycling from noncycling cells; RNA content is also a prognostic parameter in many malignancies. Two flow cytometric methods that simultaneously measure cellular RNA and DNA were developed. The 1st method is based on the use of the

metachromatic dye acridine orange, which can differentially stain these nucleic acids. This method measures total cellular or nuclear RNA content and can be used in flow cytometers that have a single source of illumination. The 2nd method uses a combination of **fluorescent dyes**, pyronin Y and Hoechst 33342. Only **double-stranded RNA** fluoresces when stained with pyronin Y. Simultaneous measurement of DNA and RNA utilizing Hoechst 33342 and pyronin Y requires the use of instruments having double excitation sources. Both methods are very sensitive to variations in dye concentration

L4 ANSWER 14 OF 15 BIOTECHNO COPYRIGHT 2005 Elsevier Science B.V. on STN
 AN 1997:27248997 BIOTECHNO
 TI Measurement of nucleic acid concentrations using the DyNA Quant(TM) and the GeneQuant(TM)
 AU Teare J.M.; Islam R.; Flanagan R.; Gallagher S.; Davies M.G.; Grabau C.
 CS J.M. Teare, Hoefer Pharmacia Biotech, 654 Minnesota Street, San Francisco, CA 94107, United States.
 E-mail: john.teare@ussfo.pharmacia.se
 SO BioTechniques, (1997), 22/6 (1170-1174), 14 reference(s)
 CODEN: BTNQDO ISSN: 0736-6205
 DT Journal; Article
 CY United States
 LA English
 SL English
 AB Molecular biology is now a routine tool in almost all biological research fields. With the exponential growth in the number of molecular biological techniques. There is a recognizable need for sensitive, accurate and precise quantitation of nucleic acids. We present here two complementary instruments designed for the quantitation of nucleic acids, the GeneQuant(TM) II and the DyNA Quant(TM) 200 Fluorometer. The GeneQuant II can rapidly determine the UV absorbance of a solution and display the calculated DNA, **RNA** or protein concentration. In addition, the GeneQuant can display calculated melting temperatures for a given DNA oligonucleotide base sequence, a useful feature for primer design. The DyNA Quant 200 quantitates DNA on the basis of the **fluorescent** Hoechst 33258 **dye/double-stranded** (ds)DNA assay. Upon binding to dsDNA, the spectral properties of the dye change such that it becomes highly fluorescent at 460 nm when excited at 365 nm. The assay has proven to be a specific and sensitive alternative method for DNA quantitation, particularly for unpurified DNA samples. Together, the GeneQuant II and the DyNA Quant 200 are a cost-effective and convenient solution to the routine protein and nucleic acid quantification needs of the molecular biologist.

L4 ANSWER 15 OF 15 BIOTECHNO COPYRIGHT 2005 Elsevier Science B.V. on STN
 AN 1993:23087012 BIOTECHNO
 TI Fluorometric assay using dimeric dyes for double- and single-stranded DNA and RNA with picogram sensitivity
 AU Rye H.S.; Dabora J.M.; Quesada M.A.; Mathies R.A.; Glazer A.N.
 CS MCB:Stanley/Donner ASU, University of California, Berkeley, CA 94720, United States.
 SO Analytical Biochemistry, (1993), 208/1 (144-150)
 CODEN: ANBCA2 ISSN: 0003-2697
 DT Journal; Article
 CY United States
 LA English
 SL English
 AB Thiazole orange homodimer (TOTO; 1,1'-(4,4,7,7-tetramethyl-4,7-diazaundecamethylene)-bis-4- ϕ 3-methyl-2, 3-dihydro-(benzo-1,3-thiazole)-2-methylidene!-quinolinium tetraiodide) and oxazole yellow homodimer (YOYO; an analogue of TOTO with a benzo-1,3-oxazole in place of the benzo-1,3-thiazole) bind with very high affinity to nucleic acids with more than a 1000-fold fluorescence enhancement upon binding. A linear dependence of fluorescence intensity on DNA concentration over a range from 0.5 to 100 ng/ml in the presence of 2×10^{-7} M TOTO or YOYO in 4 mM Tris-acetate/0.1 mM EDTA/50 mM NaCl, pH 8.2 allows sensitive quantitation of **double-stranded** DNA in a

conventional fluorometer. With nucleic acid-dye mixtures in an array of 25- μ l wells in a block of low autofluorescence plastic and detection with a laser-excited confocal fluorescence scanner, as little as 20 pg of **double-stranded** DNA can be detected per well. The array scanning method is rapid, has high throughput, and requires small amounts of sample. It also allows quantitation of single-stranded DNA and RNA.

L9 ANSWER 5 OF 7 BIOTECHNO COPYRIGHT 2005 Elsevier Science B.V. on STN
 AN 1999:29147143 BIOTECHNO
 TI Characterization of SYBR gold nucleic acid gel stain: A dye optimized for
 use with 300-nm ultraviolet transilluminators
 AU Tuma R.S.; Beaudet M.P.; Jin X.; Jones L.J.; Cheung C.-Y.; Yue S.; Singer
 V.L.
 CS V.L. Singer, Molecular Probes, Incorporated, 4849 Pitchford Avenue,
 Eugene, OR 97402, United States.
 E-mail: vicki@probes.com
 SO Analytical Biochemistry, (15 MAR 1999), 268/2 (278-288), 27 reference(s)
 CODEN: ANBCA2 ISSN: 0003-2697
 DT Journal; Article
 CY United States
 LA English
 SL English
 AB The highest sensitivity nucleic acid gel stains developed to date are
 optimally excited using short-wavelength ultraviolet or visible light.
 This is a disadvantage for laboratories equipped only with 306- or 312-nm
 UV transilluminators. We have developed a new unsymmetrical
cyanine dye that overcomes this problem. This new dye,
 SYBR Gold nucleic acid gel stain, has two fluorescence excitation maxima
 when bound to **DNA**, one centered at .sim.300 nm and one at
 .sim.495 nm. We found that when used with 300-nm transillumination and
 Polaroid black-and-white photography, SYBR Gold stain is more sensitive
 than ethidium bromide, SYBR Green I stain, and SYBR Green II stain for
 detecting **double-stranded DNA**,
 single-stranded **DNA**, and **RNA**. SYBR Gold stain's
 superior sensitivity is due to the high fluorescence quantum yield of the
 dye-nucleic acid complexes (.sim.0.7), the dye's large fluorescence
 enhancement upon **binding** to nucleic acids (.sim.1000-fold), and
 its capacity to more fully penetrate gels than do the SYBR Green gel
 stains. We found that SYBR Gold stain is as sensitive as silver staining
 for detecting **DNA**--with a single-step staining procedure.
 Finally, we found that staining nucleic acids with SYBR Gold stain does
 not interfere with subsequent molecular biology protocols.

L9 ANSWER 6 OF 7 BIOTECHNO COPYRIGHT 2005 Elsevier Science B.V. on STN
 AN 1997:27248997 BIOTECHNO
 TI Measurement of nucleic acid concentrations using the DyNA Quant(TM) and
 the GeneQuant(TM)
 AU Teare J.M.; Islam R.; Flanagan R.; Gallagher S.; Davies M.G.; Grabau C.
 CS J.M. Teare, Hoefer Pharmacia Biotech, 654 Minnesota Street, San
 Francisco, CA 94107, United States.
 E-mail: john.teare@ussfo.pharmacia.se
 SO BioTechniques, (1997), 22/6 (1170-1174), 14 reference(s)
 CODEN: BTNQDO ISSN: 0736-6205
 DT Journal; Article
 CY United States
 LA English
 SL English
 AB Molecular biology is now a routine tool in almost all biological research
 fields. With the exponential growth in the number of molecular biological
 techniques. There is a recognizable need for sensitive, accurate and
 precise quantitation of nucleic acids. We present here two complementary
 instruments designed for the quantitation of nucleic acids, the
 GeneQuant(TM) II and the DyNA Quant(TM) 200 Fluorometer. The GeneQuant II
 can rapidly determine the UV absorbance of a solution and display the
 calculated **DNA**, **RNA** or protein concentration. In
 addition, the GeneQuant can display calculated melting temperatures for a
 given **DNA** oligonucleotide base sequence, a useful feature for
 primer design. The DyNA Quant 200 quantitates **DNA** on the basis
 of the **fluorescent Hoechst 33258 dye/double**
-stranded (ds)DNA assay. Upon **binding** to
 dsDNA, the spectral properties of the dye change such that it becomes
 highly fluorescent at 460 nm when excited at 365 nm. The assay has proven
 to be a specific and sensitive alternative method for **DNA**
 quantitation, particularly for unpurified **DNA** samples.

Together, the GeneQuant II and the DyNA Quant 200 are a cost-effective and convenient solution to the routine protein and nucleic acid quantification needs of the molecular biologist.

L12 ANSWER 3 OF 4 CAPLUS COPYRIGHT 2005 ACS on STN

AN 1999:97867 CAPLUS

DN 130:322555

TI Interaction of cyanine dyes with nucleic acids. 3. The use of new cyanine dyes Cyan 13 and Cyan 40 for detection of nucleic acids in agarose gel

AU Yarmoluk, S. M.; Dubey, Igor Ya.

CS Inst. Mol. Biol. Genet., Nat. Acad. Sci. Ukr., Kiev, 252143, Ukraine

SO Biopolimery i Kletka (1997), 13(5), 419-421

CODEN: BIKLEK; ISSN: 0233-7657

PB Institut Molekulyarnoi Biologii i Genetiki NAN Ukrainy

DT Journal

LA English

AB Detection of **double-stranded** DNA (dsDNA) and single-stranded DNA (ssDNA) and **RNA** with two new **cyanine dyes** Cyan 13 and Cyan 40 is reported. Cyan 13 and Cyan 40 bind to nucleic acids to form a stable fluorescent complexes and can be used for the detection of DNA and RNA samples separated by gel electrophoresis. Sensitivity of detection is comparable to that for ethidium bromide (EtBr), a common nucleic acid staining dye.

RE.CNT 7 THERE ARE 7 CITED REFERENCES AVAILABLE FOR THIS RECORD

ALL CITATIONS AVAILABLE IN THE RE FORMAT

L12 ANSWER 4 OF 4 BIOTECHNO COPYRIGHT 2005 Elsevier Science B.V. on STN

AN 2002:34214365 BIOTECHNO

TI Stains, labels and detection strategies for nucleic acids assays

AU Kricka L.J.

CS Dr. L.J. Kricka, Department of Pathology, University of Pennsylvania, Philadelphia, PA 19104, United States.

E-mail: kricka@mail.med.upenn.edu

SO Annals of Clinical Biochemistry, (2002), 39/2 (114-129), 165 reference(s)

CODEN: ACBOBU ISSN: 0004-5632

DT Journal; General Review

CY United Kingdom

LA English

SL English

AB Selected developments and trends in stains, labels and strategies for detecting and measuring nucleic acids (DNA, **RNA**) and related molecules [e.g. oligo(deoxy)-nucleotides, nucleic acid fragments and polymerase chain reaction products] are surveyed based on the literature in the final decade of the 20th century (1991-2000). During this period, important families of **cyanine dyes** were developed for sensitive detection of **double-stranded** DNA, single-stranded DNA, and oligo(deoxy)nucleotides in gels and in solution, and families of energy transfer primers were produced for DNA sequencing applications. The continuing quest for improved labels for hybridization assays has produced a series of candidate labels including genes encoding enzymes, microparticles (e.g. quantum dots, nanocrystals, phosphors), and new examples of the fluorophore (e.g. **cyanine dyes**) and enzyme class of labels (e.g. firefly luciferase mutants). Label detection technologies for use in northern and southern blotting assays have focused on luminescent methods, particularly enhanced chemiluminescence for peroxidase labels and adamantyl 1,2-dioxetanes for alkaline phosphatase labels. Sets of labels have been selected to meet the demands of multicolour assays (e.g. four-colour sequencing and single nucleotide primer extension assays). Non-separation assay formats have emerged based on fluorescence polarization, fluorescence energy transfer (TaqMan.TM., molecular beacons) and channelling principles. Microanalytical devices (microchips), high-throughput simultaneous test arrays (microarrays, gene chips), capillary electrophoretic analysis and dipstick devices have presented new challenges and requirements for nucleic acid detection, and fluorescent methods currently dominate in many of these applications.